


SALINITY REFRACTOMETER

TECHNICAL DATA



FEATURES

- Salinity refractometers are useful for measuring the concentration of natural or artificial seawater in aquariums, holding tanks, and breeding ponds
- They can also be used for controlling the salt-water concentration when preparing pickles, measuring the concentration of saltwater used in rinsing and cooking sea products

SPECIFICATIONS

Range: 0 to 28%

Resolution: 0.2%

Accuracy: $\pm 0.2\%$

ORDERING INFORMATION

R9600 Salinity Refractometer
RPDPA1 Lens Cover
CERTIRM1 NIST Calibration Certificate