

PH050 Mini pH electrode

Sealed gel reference solution and protective bulb cover

Features:

- Rugged, polycarbonate construction
- Glass pH sensing bulb surrounded by protective teeth
- 0 to 14pH range
- 32 to 176°F (0 to 80°C) operating temperature
- Includes 39" (1m) cable with BNC connector and storage solution bottle
- 6 month warranty



Specifications

pH Range	0 to 14.0pH
Operating Temperature	32 to 176°F (0 to 80°C)
Isopotential Point	~7.0pH
Material	Polymer housing with glass bulb
Reference Junction	Ceramic
Reference Type	Ag/AgCl (Silver/Silver Chloride)
Cable Length	39" (1m)
Connector	BNC
Dimensions	5 x 0.4" (126 x 10mm)

Ordering Information:

PH050 Mini pH Electrode